

PUBLIC HEALTH SERVICE
U.S. FOOD AND DRUG ADMINISTRATION
SHELLFISH PROGRAM IMPLEMENTATION BRANCH
SHELLFISH SAFETY TEAM
5100 PAINT BRANCH PARKWAY
COLLEGE PARK, MD 20740-3835
TEL. 301-436-2151/2147 FAX 301-436-2672

SHELLFISH LABORATORY EVALUATION CHECKLIST

LABORATORY:

ADDRESS:

TELEPHONE:

FAX:

EMAIL:

DATE OF EVALUATION:

DATE OF REPORT:

LAST EVALUATION:

LABORATORY REPRESENTED BY:

TITLE:

LABORATORY EVALUATION OFFICER:

SHELLFISH SPECIALIST:

REGION:

OTHER OFFICIALS PRESENT:

TITLE:

Items which do not conform are noted by:

C – Critical K - Key O - Other NA - Not Applicable Conformity is noted by
a “√”

PART I – QUALITY ASSURANCE

Code	Item Description
	1.1 Quality Assurance (QA) Plan
<u>K</u>	<ol style="list-style-type: none"> 1. <u>Written Plan adequately covers all the following: (check \sqrt those that apply)</u> <ol style="list-style-type: none"> a. <u>Organization of the laboratory.</u> b. <u>Staff training requirements.</u> c. <u>Standard operating procedures.</u> d. <u>Internal quality control measures for equipment, calibration, maintenance, repair and performance.</u> e. <u>Laboratory safety.</u> f. <u>Internal performance assessment</u> g. <u>External performance assessment</u>
<u>C</u>	2. QA Plan is implemented.
	1.2 Work Area
<u>O</u>	1. <u>Adequate for workload and storage.</u>
<u>O</u>	2. <u>Clean and well lighted.</u>
<u>O</u>	3. <u>Adequate temperature control.</u>
<u>O</u>	4. <u>All work surfaces are nonporous and easily cleaned.</u>
	1.3 Laboratory Equipment.
<u>O</u>	1. <u>The pH meter has a standard accuracy of 0.1 unit.</u>
<u>K</u>	2. <u>pH paper in the appropriate range (i.e. 1-4) is used with minimum accuracy of 0.5 pH units.</u>
<u>K</u>	3. <u>pH electrodes consist of pH half cell and reference half cell or equivalent combination electrode (free from Ag/AgCl or contains an ion exchange barrier to prevent passage of Ag ions into the medium that may result in inaccurate pH readings).</u>
<u>K</u>	4. <u>pH meter is calibrated daily or with each use. Records maintained.</u>
<u>K</u>	5. <u>Effect of temperature has been compensated for by an ATC probe or by manual adjustment.</u>
<u>K</u>	6. <u>A minimum of two standard buffer solutions (2 & 7) are used to calibrate the pH meter. Standard buffer solutions are used once and discarded.</u>
<u>K</u>	7. <u>Electrode efficiency is determined daily or with each use following either slope or millivolt procedure.</u>
<u>K</u>	8. <u>The balance provides a sensitivity of at least 0.0001 g at a load of 5 grams.</u>
<u>K</u>	9. <u>The balance calibration is checked monthly using NIST class S, ASTM class 1 or 2 weights or equivalent. Records maintained.</u>
<u>K</u>	10. <u>Refrigerator temperature is maintained between 0 and 4°C.</u>
<u>K</u>	11. <u>Refrigerator temperature is monitored at least once daily. Records maintained.</u>
<u>K</u>	12. <u>Freezer temperature is maintained at -20°C or below.</u>
<u>O</u>	13. <u>Freezer temperature is monitored at least once daily. Records maintained.</u>
<u>O</u>	14. <u>All glassware is clean.</u>
<u>K</u>	<ol style="list-style-type: none"> 15. <u>High performance liquid chromatography system equipped with the following:</u> <ol style="list-style-type: none"> a. <u>Low dead-volume,</u> b. <u>binary solvent system delivering a pulse-free flow of 0.5-2.0 mL/min,</u> c. <u>solvent degasser,</u> d. <u>autosampler with loop suitable for 5-30 μL injections,</u> e. <u>temperature controlled column compartment capable of controlling temperature between 10 – 50°C, and</u> f. <u>fluorescence detector able to achieve the required sensitivity at excitation λ=330nm and emission λ=390nm.</u>
<u>K</u>	16. <u>Post-column reaction system equipped with the following:Reactor module</u>

	capable of maintaining 85°C, b. <u>dual reagent pumps capable of delivering accurate flows of 0.4 mL/min, and</u> c. <u>knitted reaction coil, 1 mL volume, 5 m x 0.5 mm.</u>
<u>K</u>	<u>17. Autopipettors are calibrated annually. Records maintained.</u>
<u>K</u>	<u>18. Boiling water bath with sufficient volume to cover sample/acid mixture.</u>
<u>K</u>	<u>19. Centrifuge capable of holding 50 mL polypropylene tubes and generating ~ 3000 RCF.</u>
<u>K</u>	<u>20. Microcentrifuge capable of generating ~16000 RCF.</u>
	<u>1.4 Reagents and Reference Solution Preparation and Storage</u>
<u>O</u>	<u>1. All solvents and reagents used are analytical or LC grade materials.</u>
<u>K</u>	<u>2. Water is distilled or deionized and exceeds 0.5 megaohm resistance or is less than 2 µSiemens/cm conductivity at 25°C to be tested and recorded monthly for resistance or conductivity.</u>
<u>O</u>	<u>3. Water is analyzed for residual chlorine monthly and is at a nondetectable level (<0.1 ppm) Records maintained.</u>
<u>K</u>	<u>4. Water is free from trace (< 0.5 mg/l) dissolved metals specifically, Cd, Cr, Cu, Ni, Pb, and Zn as determined annually with total heavy metal content < 1.0 mg/l. Records maintained.</u>
<u>O</u>	<u>5. Water contains < 1000 CFU/ml as determined monthly using the heterotrophic plate count method. Records maintained.</u>
<u>O</u>	<u>6. Reagents are properly stored and labeled with the date of receipt and date opened.</u>
<u>C</u>	<u>7. 0.5 M 1-heptane sulphonate is prepared the day of use or refrigerated.</u>
<u>C</u>	<u>8. pH of mobile phases and oxidant are as follows and records maintained:</u> <u>a. GTX/STX toxins mobile phase A&B is 7.1,</u> <u>b. C toxins mobile phase A is 5.8, and</u> <u>c. Oxidant is 7.8.</u>
<u>K</u>	<u>9. Mobile phases and post-column reagents are filtered through 0.2 µm nylon filter membrane before use.</u>
<u>C</u>	<u>10. Only certified reference materials are used for standard solutions. Source of the reference standard:</u>
<u>K</u>	<u>11. All primary standards are stored appropriately as per supplier recommendations.</u>
<u>K</u>	<u>12. Standards are prepared gravimetrically using “Class A” glassware.</u>
<u>K</u>	<u>13. Intermediate mixes of primary standards are made up in 0.003 M HCl (GTX/STX toxins) or Milli-Q water (C toxins), and stored appropriately.</u>
<u>K</u>	<u>14. Working standards are made up from primary standard mixes by dilution with toxin-free, deproteinated mussel or oyster extract (GTX/STX toxins) or Milli-Q water (C toxins).</u>
<u>K</u>	<u>15. Working standards are stored in the refrigerator at 4°C.</u>
	<u>1.5 Collection and Transportation of Samples</u>
<u>O</u>	<u>1. Shellstock are collected in clean, waterproof, puncture resistant containers.</u>
<u>K</u>	<u>2. Samples are appropriately labeled with the collector’s name, type of shellstock, the source, the harvest area, time, date and place (if market sample) of collection.</u>
<u>K</u>	<u>3. Immediately after collection, shellstock samples are placed in dry storage between 0 and 10°C until analyzed.</u>
<u>K</u>	<u>4. The time from collection to completion of the assay should not exceed 24 hours. However, if there are significant transportation delays, then shellstock</u>

	<p><u>samples are processed immediately as follows (circle the appropriate choice):</u>Washed, shucked, drained, frozen until extracted;</p> <p><u>b. Washed shucked, drained, homogenized and frozen;</u></p> <p><u>c. Washed, shucked drained, extracted, the supernatant decanted and refrigerated (best choice); or</u></p> <p><u>d. The laboratory has an appropriate contingency plan in place to handle samples which can't be analyzed within 24 hours due to transportation issues.</u></p>
	<p><u>5. Frozen shucked product or homogenates are allowed to thaw completely and all liquid is included as part of the sample before being processed futher.</u></p>
<p><u>PART II – EXAMINATION OF SHELLFISH FOR PSP TOXINS</u></p>	
<p><u>2.1 Preparation of Sample</u></p>	
<u>C</u>	<p><u>1. At least 12 animals are used per sample or the laboratory has an appropriate contingency plan for dealing with non-typical species of shellfish.</u></p>
<u>O</u>	<p><u>2. The outside of the shell is thoroughly cleaned with fresh water.</u></p>
<u>O</u>	<p><u>3. Shellstock are opened by cutting the adductor muscles.</u></p>
<u>O</u>	<p><u>4. The inside surfaces of the shells are rinsed with fresh water to remove sand and other foreign materials.</u></p>
<u>O</u>	<p><u>5. Shellfish meats are removed from the shell by separating the adductor muscles and tissue connecting at the hinge.</u></p>
<u>K</u>	<p><u>6. Damage to the body of the mollusk is minimized in the process of opening.</u></p>
<u>O</u>	<p><u>7. Shucked shellfish are drained on a #10 mesh sieve or equivalent without layering for 5 minutes.</u></p>
<u>K</u>	<p><u>8. Pieces of shell and drainage are discarded.</u></p>
<u>C</u>	<p><u>9. Drained meats or thawed homogenates are blended at high speed until homogenous (60-120 seconds).</u></p>
<p><u>2.2 Digestion of Sample</u></p>	
<u>K</u>	<p><u>1. Sample homogenates are extracted as soon as possible (same day) or stored in the freezer.</u></p>
<u>K</u>	<p><u>2. Sample homogenate is extracted in a 1:1 w/v ratio with 0.1 M HCl, preferably 5g tissue in 5mL acid</u></p>
<u>K</u>	<p><u>3. Homogenate/acid mixture is vortexed thoroughly before boiling to completely mix the contents.</u></p>
<u>C</u>	<p><u>4. To prevent toxin transformation, the pH of the homogenate/acid mixture before boiling is 3.0 ± 1.0, adjusted if necessary with 5M HCl or 0.1 M NaOH.</u></p>
<u>C</u>	<p><u>5. Samples are extracted in a boiling water bath for 5 minutes, in capped 50mL polypropylene centrifuge tubes.</u></p>
<u>K</u>	<p><u>6. The pH of the cooled mixture after boiling is 3.0 ± 1.0, adjusted if necessary with 5M HCl. Any sample with a pH of less than 2.0 is discarded and extracted again.</u></p>
<u>K</u>	<p><u>7. The homogenate/acid mixture is allowed to separate by gravity or by centrifugation at 2500 g for 10 minutes. Supernatant is then decanted into a scintillation vial.</u></p>
<p><u>2.3 Deproteination</u></p>	
<u>C</u>	<p><u>1. Extract is deproteinated with 30% trichloroacetic acid (50 µL TCA per 1000 µL aliquot of supernatant), vortexed thoroughly and centrifuged at 16,000 g for 5 minutes.</u></p>
<u>C</u>	<p><u>2. The pH of the deproteinated extract is adjusted to 3.0 ± 1.0 with 1.0 M</u></p>

		<u>NaOH (70 µL NaOH per 1000 µL aliquot of supernatant), vortexed thoroughly and centrifuged at 16,000 g for 5 minutes.</u>				
<u>K</u>		<u>3. An aliquot of the deproteinated, pH-adjusted supernatant is filtered through a 0.2 µm filter into two 2 mL autosampler vials (one vial for GTX/STX analysis and one vial for C-Toxins analysis).</u>				
		2.4 Assay				
<u>C</u>		<u>1. A calibration is performed upon initial instrument set up, following any major hardware maintenance activity, or when the continuing calibration verification (CCV) indicates significant drift (> 30% for individual toxin) from the calibration. Records maintained.</u>				
<u>K</u>		<u>2. For GTX/STX toxins, no more than ten samples should be made between standard analyses. For C toxins, no more than five samples injections should be made between standard analyses.</u>				
<u>K</u>		<u>3. 10 µL is injected for GTX/STX toxins and 5 µL is analyzed for C-toxins.</u>				
<u>K</u>		<u>4. Samples are stored in the sample compartment at 4°C during analysis.</u>				
<u>O</u>		<u>5. A column heater is used in the analysis.</u>				
<u>O</u>		<u>6. The appropriate analytical column is used.</u> <u>a. GTX/STX Toxins: Zorbax Bonus-RP column, 4.6 mm x 150 mm, 3.5 µm, Agilent catalog number 863668-901 or equivalent.</u> <u>b. C Toxins: BetaBasic 8, 4.6 mm x 250 mm, 5 µm, Fisher catalog number 71405-254630 or equivalent.</u>				
		2.5 System Suitability				
<u>K</u>		<u>1. The correlation coefficient for the linear regression (r²) must be > 0.990 for each individual toxin.</u>				
<u>K</u>		<u>2. Resolution and Retention Time Criteria.</u> <u>a. GTX/STX Toxins.</u> <u>i. Matrix peak must be at least 70% baseline resolved between GTX3 and GTX2.</u> <u>ii. GTX5 must be at least 40% baseline resolved between dcGTX3 and dcGTX2.</u> <u>iii. dcSTX and STX must be at least 70% baseline resolved.</u> <u>iv. GTX4 retention time should be between 5 and 7 minutes.</u> <u>b. C Toxins.</u> <u>i. C1 and C2 must be at least 70% baseline resolved.</u> <u>ii. C1 retention time should be between 5 and 8 minutes.</u>				
		2.6 Calculation of Toxicity				
<u>C</u>		<u>1. The toxicity of the individual toxins is calculated as follows:</u> $\mu\text{gSTXdiHCle q}/100\text{g} = \mu\text{M} \times \frac{372.2}{1000\text{mL}} \times \frac{\text{Fvol}}{\text{Ext.vol}} \times \left(\frac{\text{Wt} + \text{Vol}}{\text{Wt}} \right) \times \text{ReTx} \times 100$ <hr/> <u>Where:</u> <u>µM = Concentration of toxin in the extract, in µM;</u> <u>Fvol = Final volume of the deproteinized extract (1120 µL);</u> <u>Ext.vol = Volume of crude extract used (1000 µL);</u> <u>Wt = Weight of sample used;</u> <u>Vol = Volume of acid extractant used (e.g. 5 mL); and</u> <u>ReTx = Relative toxicity of toxin vs. Saxitoxin.</u>				
		<u>Relative Toxicity Values</u>				
		<table border="1"> <tr> <td><u>Toxin</u></td> <td><u>ReTx</u></td> <td><u>Toxin</u></td> <td><u>ReTx</u></td> </tr> </table>	<u>Toxin</u>	<u>ReTx</u>	<u>Toxin</u>	<u>ReTx</u>
<u>Toxin</u>	<u>ReTx</u>	<u>Toxin</u>	<u>ReTx</u>			

		<u>GTX1</u>	<u>0.9940</u>	<u>NEO</u>	<u>0.9243</u>
		<u>GTX2</u>	<u>0.3592</u>	<u>STX</u>	<u>1.0000</u>
		<u>GTX3</u>	<u>0.6379</u>	<u>dcSTX</u>	<u>0.5131</u>
		<u>GTX4</u>	<u>0.7261</u>	<u>C1</u>	<u>0.0060</u>
		<u>GTX5</u>	<u>0.0644</u>	<u>C2</u>	<u>0.0963</u>
		<u>dcGTX2</u>	<u>0.1538</u>	<u>C3</u>	<u>0.0133</u>
		<u>dcGTX3</u>	<u>0.3766</u>	<u>C4</u>	<u>0.0576</u>

C 2. The individual toxicities for each toxin are summed to obtain the overall sample toxicity in µg STX equivalents/100 g (µg/100 g)

C 3. Any value greater than 80 µg STX equivalents /100 g of meat is actionable.

REFERENCES

1. AOAC Official Methods of Analysis (2011). AOAC Official Method 2011.02 Paralytic Shellfish Toxins in Mussels, Clams, Oysters, and Scallops Post-Column Oxidation (PCOX) Method.
2. Adams, W.N. and S.A. Furfari. 1984. Evaluation of laboratory performance of the AOAC method for PSP toxin in shellfish. *J. Assoc. Off. Anal. Chem.* Vol 67, 6:1147-1148.
3. American Public Health Association. 1970. *Recommended Procedures for the Examination of Sea Water and Shellfish*, 4th Edition. APHA, Washington, D.C.
4. American Public Health Association. 192. *Standard Methods for the Examination of Dairy Products*, 16th Edition. APHA, Washington, D.C.
5. Association of Official Analytical Chemists International. 1990. *Methods of Analysis*, 15th Edition. AOAC, Arlington, VA.
6. APHA/WEF/AWWA. 1992. *Standard Methods for the Examination of Water and Wastewater*, 18th Edition. APHA, Washington, D.C.
7. Title 21, Code of Federal Regulations, Part 58, *Good Laboratory Practice for Nonclinical Laboratory Study*. U.S. Government Printing, Washington, D.C.
8. National Research Council. 1996. *Guide for the Care and Use of Laboratory Animals*. National Academy Press, Washington, D.C.
9. Personal communication with USFDA Washington Seafood Laboratory Branch, Office of Seafood, CFSAN, 1998-1999.

<u>LABORATORY STATUS</u>	
<u>LABORATORY</u>	<u>DATE</u>
<u>LABORATORY REPRESENTATIVE:</u>	
<u>PARALYTIC SHELLFISH POISON COMPONENT: PARTS I AND II</u>	
<u>A. Results</u> <u>Total # of Critical (C) Nonconformities</u> <u>Total # of Key (K) Nonconformities</u> <u>Total # of Critical, Key, and Other (O) Nonconformities</u>	<hr/> <hr/> <hr/>
<u>B. Criteria for Determining Laboratory Status of the PSP Component</u> <p><u>1. Does Not Conform Status</u> The PSP component of this laboratory is not in conformity with NSSP requirements if:</p> <ul style="list-style-type: none"> a. <u>The total # of Critical nonconformities is ≥ 3 or</u> b. <u>The total # of Key nonconformities is ≥ 6 or</u> c. <u>The total # of Critical, Key, or Other is > 10</u> <p><u>2. Provisionally Conforms Status</u> The PSP component of this laboratory is determined to be provisionally conforming to NSSP requirements if the number of critical nonconformities is <u>≥ 1 by < 3.</u></p>	
<u>C. Laboratory Status (circle appropriate)</u> <p style="text-align: center;"><u>Does Not Conform – Provisionally Conforms – Conforms</u></p>	
<u>Acknowledgement by Laboratory Director/Supervisor:</u> <u>All corrective Action will be implemented and verifying substantiating documentation received by the Laboratory Evaluation Officer on or before</u> <hr/>	
<u>Laboratory Signature:</u> _____	<u>Date:</u> _____
<u>LEO Signature:</u> _____	<u>Date:</u> _____