

IV. Summary Table for QPCR Methods

Validation Criteria	Test	Number of Tests	Statistical Test – Data Handling	Recommended Acceptable Criteria (if available)
1. Accuracy/Trueness	1. & 2. For each shellfish of interest- Spike with test analyte: Spiked sample Blank samples Range: Low, but detectable to 10 ⁷ Run by QPCR & by plate	1. & 2. 20 samples 2 aliquots per sample: 1 spiked 1 blank Run by both plate and QPCR	1. Log transformation Avg QPCR MPN/Avg plate count X 100	
2. Measurement Uncertainty			2. Log transformation Find difference between plate and corresponding QPCR Find 95% confidence interval for differences	
3. Precision	3. & 4. For each shellfish of interest: Ten samples: Spike with three concentrations of analyte (low, medium & high). Analyze by plate method and QPCR	3. & 4. 10 Samples: 4 aliquots per sample Low, medium, high concentration + blank 2 replications of each concentration QPCR, 1 replication at each concentration plate method, 1 blank Run by plate and QPCR	3. Log transformation Nested ANOVA If F α 0.05 is not significant, compare standard deviation of total over all concentrations to standard deviation of 3 tube decimal dilution (0.317) with one sided t test α 0.05.	3. Calculated F < F α 0.05 S _L < 0.317
4. Recovery			4. Log transformation Average QPCR replicates Subtract Avg QPCR from associated plate count Perform single ANOVA	4. Calculated F < F α 0.05 If significant , use Tukey’s HSD to compare recovery by concentration

<p>5. Specificity</p>	<p>5. For each shellfish of interest: Spike with test organism (at low level)& interfering organism at moderate level Run by QPCR</p>	<p>5. 1 samples per interfering organism(io): 3 aliquots 1 spiked w analyte organism 1 spiked w analyte + io 1 blank Run: 5 replicates of spiked samples 1 repliccate of accompanying blank</p>	<p>5. Log transformation Specificity index (SI): $SI = \frac{\text{Analyte Spike}}{\text{Analyte io Spike}}$ For SI > or < 1: Perform two sided t test</p>	<p>5. $[SI_{avg} - m_{0=1}] > u$</p>
<p>6. Linear Range</p>	<p>6, 7, & 8. For each shellfish of interest: Each sample divide into 10 aliquots Spike with 9 concentrations $10^0, 10^1 \dots 10^8$ Analyze by plate count and QPCR</p>	<p>6, 7, & 8. 5 Samples: 9 spike concentrations per sample: 2 replicates by QPCR 1 replicate by plate method</p>	<p>6. Plot Critical Threshold(Ct) verses Plate Counts (log transformed) Provide equation for the line Determine if linear relationship exists , compute correlation coefficient (Pearson's r)</p>	<p>Test r for significance</p>
<p>7. Limit of Detection</p>			<p>7. If linear, use $y=mx+b$ Where: $y = \# \text{ cycles of PCR amplification}$ Using slope and y-intercept from 6. above, set $y =$ to total number of PCR cycles and solve for x Take antilog of X</p>	<p>$X = 1$ If different from 1, determine if significantly different at 95% confidence interval.</p>
<p>8. Limit of Quantitation/ Sensitivity</p>			<p>8. If $x = 1$, LOQ based on 3 tube, 3 dilution ratio of 0.01, 0.001, 0.0001</p>	

<p>9. Ruggedness</p>	<p>9. For each shellfish of interest: Two aliquots Spike with analyte –range of concentrations Prepare two separate batches(or lots) of reagents Analyze by QPCR , 1 aliquot use set 1 reagents, second aliquot use set 2 reagents</p>	<p>9. 10 Samples Use set 1 and set 2 prepared reagents for all ten samples Samples should cover range of concentrations. Process samples over a period of days</p>	<p>9. Log transformation Perform two sided t-test at 95% Confidence interval</p>	<p>9. If $\bar{X}_{A\text{ Avg}} - \bar{X}_{B\text{ Avg}} < u$ Method sufficiently rugged</p>
--	--	---	---	--